Manual de uso / User Manual

Powerlift Up 2

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Declara que el producto / *This is to certify that the product*

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Complies with the requirements of RD 1591/2009 of 16th October, that regulates medical devices and implementing Directive 93/42/ECC, modified by Directive 2007/47/EC.

Declaration in accordance with Annex VII of RD 1591/2009. Medical devices class I

Reference standard: **UNE-EN ISO 10535 "Cranes for moving people with disabilities. Requirements and test methods"**, published in October 2007.

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
1.- FEATURES, WARRANTY, TECHNICAL SERVICES

This product is intended by the manufacturer to be used on human beings for the purpose of treatment or alleviation of disease, and treatment, alleviation of, or compensation for, an injury or disability, in accordance with current Spanish regulations on medical devices.

This product may only be used for lifting, lowering and moving patients. The patient should have some mobility and capacity of fastening because the hoist is not destined for long transfer, if not to incorporate/descends patients from chair/toilet/bed for changing nappies and/or short transfer.

Our products have a 2-year warranty against any manufacturing fault (12 months in the case of batteries and electrical equipment). If the buyer or the end user allows unauthorized personnel to operate the product or fails to observe the cleaning and maintenance recommendations contained in this manual, all warranty rights and liability claims will be forfeited. For reasons of hygiene, used slings may not be returned.

Tecnimoem sells all its products through a network of distributors. In the event of a problem, the end user should contact the DISTRIBUTOR from whom they have purchased the product. Tecnimoem provides distributors with the training and spare parts needed to undertake maintenance work on the products.

TECNIMOEM 97, S.L. operates a policy of continuous improvement of all its products. Hence, device specifications may change without prior notice.

2.- SAFETY PRECAUTIONS

Any person using or assembling these items must first read and understand the assembly, operating and cleaning instructions provided.

Failure to follow these instructions may result in damage to property and/or personal injury.

2.1.- Recomendations

 **HOIST NEVER CAN BE USE BY THE OWN PATIENT**, must be used by a person with knowledge about it. Familiarize yourself with the safety controls and devices before operating the hoist. Tecnimoem

recommends using the sling supplied with the hoist, as it has been previously tested. If you use a different sling, check that it is EC marked and that its size and other characteristics are appropriate to the patient. Do not use a sling that is frayed or worn, and always follow the instructions provided with it.

Check that the combined weight of the patient and the sling does not exceed the maximum working load displayed on the hoist. (pages 6-7)

A loaded hoist should always be used on flat, smooth surfaces that are free of obstacles. If moving a patient along a sloping surface, the gradient must not be greater than 5° and it is recommended to use the assistance of a helper. Never leave a hoist loaded with a patient on a sloping surface. Transfer the patient with the hoist legs in parallel position at a maximum speed of 3 km/h (0.8 m/s).

Keep the hoist away from water, humidity, splashes and corrosive atmospheres (indoor swimming pools, the sea, etc.). Do not use the electric hoist in the shower or recharge the battery in the bathroom.

Do not force the safety controls and devices. Do not push the mast, the boom or the patient to move the hoist. Lift or lower the patient with the rear wheel brakes on. Apply the brakes only to transfer a patient from the hoist to a chair or a bed, or the reverse.

Is not allowed use the hoist without shoes.

If hoist have been displaying in the sun during long time, could have parts with risk of burn.

BEFORE MOVE THE PATIENT Consider individual hazards, including:

- Mobility and capacity of fastening of patient
- Risk of crushing/squeezing.
- Potential falls.
- Disturbed patients.
- Patients with learning difficulties.
- Small children and pregnant women.
- Persons without the mental capacity to recognize unsafe behavior.
- Unauthorized persons.

Anybody who is authorized to use the hoist must be able to do so in a safe and controlled way.

Due to the continuous, smooth lifting action of the moving parts, there is risk of entrapment. When operating these parts, the user and/or the carer must make sure that there are no body parts in the areas where they may become trapped.

2.2.- Electrical safety precautions

All the hoist's electrical connections must comply with the International Electrotechnical Commission (IEC) standards. The power supply should be equipped with an earth-leakage circuit breaker with a maximum operating current of 30 mA, as prescribed by standard IEC 364-5-53.

This product complies with the standards EN 60601-1 and EN 60601-1-2 on electrical apparatus and the electromagnetic interference of medical devices. Therefore, it does not interfere or is not interfered when combined with other medical devices that also comply with the electromagnetic standards.

Some apparatus, particularly older devices, which do not comply with electromagnetic compatibility standards may, however, cause or be affected by interference when used with this hoist. If using such apparatus, you must make sure that any potential malfunction will not harm the patient or any other person.

 Before moving the hoist, make sure the power supply cable is unplugged from the mains. Any work on electrical parts may only be undertaken by authorized, qualified personnel.

3.- ASSEMBLY

 No tools are necessary to assemble the hoist, although some parts are considerably heavy. It is therefore advisable to use the assistance of a helper when handling such parts. Put particular attention to avoid being knocked or trapped.

3.1.- Components *(see page 3)*

BASE STRUCTURE

- 1.- Base
- 2.- Back castors Ø 100mm(with brake)
- 3- Side castors Ø 50mm w/o brake
- 4.- Front castors Ø 50mm w/o brake
- 5.- Platform feet

MAST – HANDLE STRUCTURE

- 6.- Mast
- 7.- Handle
- 8.- Handle for user

- 9.- Handle for patient
- 10.- Support for sling
- 11.- Electric motor
- 12.- Battery charging base
- 13.- Removable battery
- 14- Remote control
- 15.- Locking knob (top)
- 16- Locking knob (lower)
- 17.- Mast bracket (actuator lower support)
- 18.- Boom bracket (actuator top support)
- 19.- Support for knees
- 20- Mast for knee's support

3.2.- Joining the pre-assembled structures *(see page 4)*

Remove packing carefully.

 **DON'T REMOVE THE GREY PLASTIC PIECES (A) OF THE MAST BASE**

Place the base structure on the ground with the back casters braked. Remove the top locking knob (15) and the lower locking knob (16) of mast.

Insert mast (6) into the base (1), following insert the lower locking knob (15) fixing the mast (6) to the base (1) and mast for knee's support (20).

Finally, hang lower locking knob (16) that adjusts the base (1) to the mast (6).

Knee's support (19) have to put with the mast for it (20), and you will have five different positions.

Finally, we add the platform feet

The electric actuator is not supplied and assembled on the hoist in the factory. If you need to assemble the actuator, loosen the ring (18) and the mast (17). To disassemble do the opposite operation.

Connect **FIRMLY** in the base charge the remote control and the linear actuator .

4.- OPERATING INSTRUCTIONS

4.1.- Operating

To spread the legs, step on the left pedal. To narrow the legs, step on the right pedal. While transferring patients and to negotiate doors and confined corridors, the legs should preferably be narrowed.

Lifting and lowering the patient is achieved by the electric actuator's action on the hoist boom. These movements are easy to control using the two buttons (up and down) provided on the hand control. When the boom reaches the upper or lower limit of its run, it stops automatically.

When lifting or lowering a patient, it is recommended to press and hold the relevant button. Repeatedly pressing and releasing the buttons will cause the actuator motor to start and stop each time, using more battery and shortening its operating lifetime.

The electrical system includes a **emergency stop red button**. When the red button is in the out position, the battery will operate normally. When the red button is depressed, the power from the battery is cut off. **To reset the button to the out position, turn it clockwise and release.**

The **knee's support** is adjustable: allowing five different heights from floor.

This support is basculable to guarantee the best support for knees adapting to each patient. It has a strap with safety lock in case it is necessary to adjust the legs of the patient

4.2.- Recharging battery

System has a noisy system when the battery is near to finish, with this is possible made complete cycles. If you only have one removable battery and one charging base, place the hoist near an electric socket, connect the power supply cable to the charging base and plug into the mains (110-220V

While the battery is charging the actuator fails, the emergency stop should not be pressed and green LEDs "on" lights. The battery indicator lines flash. When fully charged the 4 battery indicator lines will remain fixed The green light "on" will remain on until you disconnect from the net-work.

If the hoist is used daily, it is recommended to recharge the battery during the night. The charger will not allow the battery to overcharge.

 When not using the hoist, it is recommended to re-charge the battery at least once every three months to ensure it is always in good working condition.

 Do not disconnect the black charger jack plug by pulling the cable. When not recharging the battery, it is recommended to keep the mains cable plugged into the charging base to prevent any liquids accidentally entering the socket.

If you have an **extra charging base** (on request), you can remove the battery from the charging base mounted on the hoist and charge it using the extra charging base, suitably located near a mains electricity socket.

4.3.- Changing the battery

To remove the battery from the charging base, hold the upper handle, press and remove carefully of its support. To place the battery back in the charging base, slide it on position, join the upper handle on the support until you hear a click

4.4.- Overload

Control box beeps twice if the actuator is overloaded (above the current limit) Start work after some minutes when the overload will be eliminated. Work cycle aprox: 2 minutes in continuous use and 18 minutes in break.

4.5.- Improper use

Improper use of this hoist may result in injury to the users and/or damage to the device, in addition to invalidating the warranty. Examples of improper use are.

- Use by persons who have not read this user manual and/or who have not been trained by a member of the authorized personnel.
- Use of the hoist, functions, accessories or movement by persons not qualified to operate the device safely.
- Use of the electrical functions by more than one person simultaneously.
- Use with a load exceeding the maximum working load displayed on the hoist.
- Faulty connection to the mains when charging the battery.
- Connecting any devices to the hoist that have not been authorized by the manufacturer.
- Pulling the power supply cable to move the hoist.
- Cleaning the hoist with excessive water, pressure jets or in a wash-tunnel.
- Using the hoist outdoors or to move a patient inside a vehicle.
- Using the hoist on soft, unprepared terrain.
- Using the hoist on terrain with a slope gradient greater than 5° (loaded with a patient).
- Extreme, intensive use of the actuator failing to observe the maximum working load displayed on the label.
- Using devices and/or accessories other than those recommended by the manufacturer.
- Any other use of the hoist that is not in accordance with its intended purpose.
- Do not respect work cycles: 2 minutes in continuous use and 18 minutes in rest.

5.- SLING

Tecnimoem recommends using the sling supplied with the hoist, which may be used with patients having a high degree of physical dependency provided that they have good head control and some mobility in inferior part. This sling allows easy access by technical and assistant personnel and may be used both on a sitting and a recumbent patient.

The sling is 100% polyester. The materials used in its manufacture are hypoallergenic. No skin reactions have been reported, but if any should appear it is recommended to consult with the prescriber.

Three different heights with security support.

5.1.- Sling safety and maintenance

Choosing the right size of sling and correctly fitting it on the area requiring treatment are crucial to ensure its purpose is fulfilled safely and efficiently, with a satisfactory lifetime. Bear in mind that a loose sling may diminish its therapeutic efficacy and an excessively tight sling may have an adverse effect on blood circulation and/or illnesses other than those you set out to treat.

SAFETY Never use with loads exceeding the maximum working load of the hoist (as displayed on the hoist)

Before using a sling on a patient, check that it is in perfect working condition. It is recommended to use it previously on healthy persons to become familiar with its use.

Do not use the sling if it appears torn, worn, frayed or distorted.

Do not use the end blue loops – these are safety attachments in case any of the other loops should break.

MAINTENANCE Recommended 30°C wash using neutral soap. Do not bleach or dry-clean. Do not iron. Dry on a flat surface. Do not expose to direct sunlight or any source of heat. Do not tumble dry.

5.2.- Positioning of sling to patient, sitting in a chair, wheelchair or lateral of bed.

Inclining the patient to the front and checking if he is well-fastened to reduce risk of fall and to be the patient secure.

Put sling under the back of patient, in the middle, with the opening in symmetrical position with the axis of vertebral column.

Adjust the belt for better security.

5.3.- Fastening of sling to patient handle

Move the hoist closer with opened legs around the patient, in such a way he can put his feet in the platform without risk of injury. If would be necessary, adjust the belt of the knees.

With the caster blocked, put off arms of patient outside the sling's tapes and fastening this into the handle, and choose the correct height (red, green or yellow).

For make it easier, is better if the arm of hoist is in his low position avoiding strain in the tapes.

5.4.- Raising up of patient

 Before lifting the patient, recliner patient's back and check if the sling is correctly fastened. If some element is not correct, take down the patient and correct the problem.

 Is recommendable blocked back caster when you are lifting a patient, to reduce risk of knock over the hoist.

Use the remote control for lifting the patient.

5.5.- Transport of patient

 During the transport, with the patient in sling, DO NOT slide the hoist in irregular floor because could knock over the hoist.

ALWAYS use the handle to move the hoist. Move patient with the legs closed.

5.6.- Descent patient

Move hoist until (wheelchair with their brakes activated, chair, toilet, bed...) where you want put the patient.

Use remote control for lifting.

 Be sure that the patient is in secure place.

6.- CLEANING AND MAINTENANCE

6.1.- Cleaning and disinfecting

- Make sure the hoist is immobilized and disconnected from the electricity mains.
- Press the red stop button to shut off the electrical functions.
- Never clean the hoist with copious water, with a high-pressure jet or in a wash-tunnel.
- Do not use water at a temperature above 60° C.

- Avoid any excess water on the connector sockets.
- Dry the hoist thoroughly before using again.
- Use clean cloths, disposable or recyclable. .
- Use a diluted solution of detergent or disinfectant, or a spray-on disinfectant.

6.2.- Safety rules for maintenance

Maintenance operations may only be performed by duly qualified and authorized technical personnel. Contact your local Tecmoem 97, S.L. distributor for any specific issue in connection with maintenance and/or servicing.

Before performing any maintenance or repair operation:

- Make sure the hoist is immobilized (if no movement is envisaged) and disconnected from the electricity mains.
- Press the red stop button to shut off the electrical functions.
- Do not under any circumstances open or puncture an electric motor.

6.3.- Preventative Maintenance

It is advisable to carry out a yearly maintenance check of the hoist and its accessories to ensure they are kept in good working order. Special attention should be paid to the following aspects:

- Operating controls and mechanisms.
- The hinges where the hoist and accessories move.
- The condition of the electric cables and the water-tightness of the electrical devices.
- The condition of the sling (wear and fraying, tears, distortion, etc.).

It is recommended to grease the hoist's hinges regularly.

Adapt the frequency and intervals of maintenance checks to the condition of the hoist and the specific circumstances of its use.

6.4.- List of spare parts

On request, electrical equipment (battery, actuator, control unit, charging cable or control box), harness, wheels and other components of the metal structure of the crane can be replaced.

6.5.- Transport and storage

During transport, the hoist must be in the low position, with its functions disconnected and the brakes on the wheels on. Additionally, it must be secured with straps and protected against water and humidity (75%), at a temperature between -20°C and +50°C.

7.- TROUBLESHOOTING

PROBLEM	CAUSE	REMEDY
The actuator is not working	The red emergency stop is activated	Turn the red button clockwise to reactivate
	The battery is connected to the electricity mains for recharging	Unplug the black power supply
	The battery is low	Recharge the battery
	The battery is not properly connected	Remove the battery and replace it
	Control or motor cables wrong connected	Insert the cables correctly to the bottom
	Hand control and/or actuator failure	Replace the hand control and/or the actuator
The battery will not charge	Mains failure	Check mains supply
	Malfunctioning battery and/or charging base	Replace the battery and/or the charging base
The actuator stops	Hoist overloaded	Reduce the load and wait a few minutes
	The battery is low	Recharge the battery
Noisy hinges	Lubrication is requires	Apply lubricant to the hinges
The wheels do not rotate or brake properly	Dirty	Clean the wheels

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CE

TECNIMOEM 97, S.L.
Ctra. Logroño - Mendavia, Km. 5
31230 Viana - Navarra – España (Spain)
T.: (+34) 948 646 213 - FAX: (+34) 948 646 380
www.tecnimoem.com